

# Feline Haptoglobin (Hp) TurboReader™ Assay



**DETECTING  
SUBCLINICAL DISEASE  
AND MONITORING  
LONG-TERM  
INFLAMMATION  
ACTIVITY  
IN CATS**



## **Feline Hp TurboReader™ Assay**

Sample type: ..... Serum or Li-heparin Plasma  
Sample volume: ..... 20 µl  
Measuring range: ..... 17 – 1000 mg/dl  
Kit size: ..... 20 tests  
Analysis time: ..... 3 minutes  
Clinical decision limit: ..... 300 mg/dl

Article No. 2532-01

Feline specific calibrated assay.

**When SAA does not agree with the clinical picture, follow-up with Hp.**